

Intralock


Titanium Femur



TRAUFIX
READY TO RUN

INTRALOCK LATERAL FEMORAL NAIL


MULTIPLE LOCKING OPTIONS



Innovative distal locking device for precise collocation of 4.5mm screw



TO PERFORM LOCKING WITH DEVICE
PRELIMINARY STEP. Calibration of distal strip to generate hole of distal locking device




MEASURES

Intralock Lateral Femoral Nail

Length	Diameter 9mm	Diameter 10mm	Diameter 11mm
200mm	Code 157.20	Code 158.20	Code 216.20
240mm	Code 157.24	Code 158.24	Code 216.24
320mm	Left Code 136.32 / Right Code 138.32	Left Code 137.32 / Right Code 139.32	Left Code 180.32 / Right Code 181.32
340mm	Left Code 136.34 / Right Code 138.34	Left Code 137.34 / Right Code 139.34	Left Code 180.34 / Right Code 181.34
360mm	Left Code 136.36 / Right Code 138.36	Left Code 137.36 / Right Code 139.36	Left Code 180.36 / Right Code 181.36
380mm	Left Code 136.38 / Right Code 138.38	Left Code 137.38 / Right Code 139.38	Left Code 180.38 / Right Code 181.38
400mm	Left Code 136.40 / Right Code 138.40	Left Code 137.40 / Right Code 139.40	Left Code 180.40 / Right Code 181.40

4.5mm LOCKING SCREW FOR INTRAMEDULLARY NAIL INTRALOCK (CODE 166, LENGTH 25mm-80mm)
6.5mm SUBTROCHANTERIC LAG SCREW INTRALOCK (CODE 167, LENGTH 60mm-105mm)
END CAP FOR FEMORAL NAIL INTRALOCK (CODE 168.12)

INTRALOCK RETROGRADE FEMORAL



Code 166



Code 179



Code 178



Code 168.12

MEASURES

Intralock Retrograde Femoral

Length	Diameter 9mm	Diameter 10mm	Diameter 11mm
240mm	Code 161.24	Code 159.24	Code 160.24
260mm	Code 161.26	Code 159.26	Code 160.26
320mm	Code 176.32	Code 174.32	Code 175.32
340mm	Code 176.34	Code 174.34	Code 175.34
360mm	Code 176.36	Code 174.36	Code 175.36
380mm	Code 176.38	Code 174.38	Code 175.38

- 4.5mm LOCKING SCREW FOR INTRAMEDULLARY NAIL INTRALOCK (CODE 166, LENGTH 25mm-80mm)
- 6.4mm SCREW FOR RETROGRADE NAIL INTRALOCK (CODE 179, LENGTH 25mm-120mm)
- 6.4mm COMPRESSION SCREW FOR RETROGRADE NAIL INTRALOCK (CODE 178, LENGTH 40mm-120mm)
- WASHER OF 4.5mm - 6.5mm (CODE 182.65)
- END CAP FOR FEMORAL NAIL INTRALOCK (CODE 168.12)

Titanium Alloy 6AL4V ELI ASTM F136

Special for surgical implants

Design by **TRAUFIX**
READY TO RUN




MADE IN MEXICO



FEMORAL NAIL SET

INSTRUMENTATION AND SET FOR ANTEROGRADE AND RETROGRADE NAILS. ALL-IN-ONE

QTY.	EQUIPMENT
1	Initiator punch
1	Universal puller
1	Allen key
1	Large allen key
1	Quick release T-handle
1	Guide Wire Gauge 120
1	Graduated drill with stop \varnothing 6.4
1	Depth gauge 115 mm
6	Guide Wire (3 Threaded / 3 Smooth Tips) 2.5 x 35cm
2	Nail Set Screw
1	U-lock strip
1	Guide Wire Positioner
1	Block for distal probe
2	Drill sleeve exterior \varnothing 8.0
1	Screwdriver for sliding screw \varnothing 6.5
1	Drill sleeve short \varnothing 8.0
1	Quick Release Hex Driver Bit for Nail clamping \varnothing 6.5
1	4.5 Bolt T-Handle Hex Driver
1	Stylus with T-handle
1	Short Tissue Dilator (Awl) \varnothing 8.0
1	T-handle drill bit \varnothing 5.0
1	Long Tissue Dilator (Awl) \varnothing 8.0
1	Drill sleeve interior \varnothing 2.5
1	Drill sleeve interior \varnothing 6.4
1	Drill sleeve interior \varnothing 4.3
1	Drill sleeve interior \varnothing 4.3
1	Drill sleeve interior \varnothing 5.2
1	Drill with stop 4.3 x 30cm
1	Drill Bit 4.3 x 30cm
1	T-Handle Hex Driver
1	Drill 5.2 x 25 cm
10	Rhying Tip Set \varnothing 8.5 \varnothing 9.0 \varnothing 9.5 \varnothing 10.0 \varnothing 10.5 \varnothing 11.0 \varnothing 11.5 \varnothing 12.0 \varnothing 12.5 \varnothing 13.0
1	Impactor / striker
3	Screw for strips
1	Probe / Distal Locking Devices
1	Fabric protector
1	Tramodular guide with T-handle
1	Flexible Rhyme
1	Starter Drill with Stop \varnothing 13
1	Distal strip 180-240
1	Strip for proximal and subtrochanteric locking (Block to femoral head)
1	Distal strip 320-440
1	2.5 Olive Tip Guide Wire
1	Drill 6.3
1	Shirt 6.3
1	Guide without Olive
1	Addition
1	Screw for attachment
1	Special strip
1	Retrograde strip
1	Satelite



TRAUFIX

internationalsales@traufix.com.mx
(+52) 55 8106 3240

www.traufix.com

DIPROMEDIC MEXICO

customerservice@dipromedic.com
(55) 9000 4888/ 89/ 90/ 91

DIPROMEDIC PERU

gustavo.dominguez@dipromedic.com
ventas.interdipro.peru@dipromedic.com
(+51) 947 732 069

 traufix.mx

 traufix